

## Titration accuracy skills check list

Tick off each criterion on the check list as you carry out your titration

Rinse ALL glassware with distilled water	
Rinse burette and pipette with solutions they will contain	
Read burette and pipette from meniscus, at eye level	
No air bubbles present in pipetted aliquot	
Swirl conical flask during titres	
Use distilled water to rinse off residue on the sides of the flask	
Obtain a very faint pink colour at the equivalence point	
Record 3 concordant titres which all agree within 0.2mL	

## Titration accuracy skills check list

Tick off each criterion on the check list as you carry out your titration

Rinse ALL glassware with distilled water	
Rinse burette and pipette with solutions they will contain	
Read burette and pipette from meniscus, at eye level	
No air bubbles present in pipetted aliquot	
Swirl conical flask during titres	
Use distilled water to rinse off residue on the sides of the flask	
Obtain a very faint pink colour at the equivalence point	
Record 3 concordant titres which all agree within 0.2mL	